

Using InDraw software to Rapidly Convert Image Structures to SMILES

1. Software Installation

Download and install the InDraw software via the <https://www.integle.com/static/indraw>.



2. Structure Localization

Upload the PDF document and use PSprompt to locate the target chemical structure in large language models such as ChatGPT.

2025_Dyad System of B... PDF

2025_Dyad System of B... PDF

You are a **photosensitizer** literature information **extraction assistant**. Your task is to **read** all uploaded literature, including the main text and supporting information, and **extract** the name, type, and structural location of each photosensitizer molecule mentioned. You must output the data in **Markdown Table** format. The table must contain 4 columns. If any required information is not mentioned in the text, fill it with "N/A". Only use information from the provided document. Do NOT use external knowledge or make assumptions.

Classification Rules:
Type I : if ANY of the following conditions are met:
Explicit generation of: $O_2^{\bullet-}$, $\cdot OH$, or radicals
ROS detected using DHR123, HPF, EPR, or similar radical-detecting probes

Type II :
Explicit generation of 1O_2
 1O_2 is detected using ABDA, SOSG, or other singlet oxygen probes
No evidence of $O_2^{\bullet-}$, $\cdot OH$, or other radicals is reported

Uncertain: If ROS type cannot be clearly determined from the text, mark as N/A

The Markdown table must contain exactly 4 columns in this order: | PS name | Type (I or II) | Type of ROS | Location |

Column Definitions
PS name: Name or abbreviation exactly as written in the text
Type (I or II): type I / type II / N/A
Type of ROS: $O_2^{\bullet-}$ / $\cdot OH$ / 1O_2 / $O_2^{\bullet-}$, $\cdot OH$, 1O_2 / N/A
Location: Position of the photosensitizer structure in the document (e.g., Figure 2)

If a photosensitizer produces both Type I and Type II ROS, it must still be classified as Type I .

Human


Based on the provided main article and supporting information, the following photosensitizers (PSs) were identified:

PS name	Type (I or II)	Type of ROS	Location
4a	I	$O_2^{\bullet-}$, 1O_2	Figure 1a (main text); Figure 2 (ROS generation)
4b	I	$O_2^{\bullet-}$, 1O_2	Figure 1a (main text); Figure 2 (ROS generation)
4c	I	$O_2^{\bullet-}$, 1O_2	Figure 1a (main text); Figure 2 (ROS generation)
5a	I	$O_2^{\bullet-}$	Figure 1a (main text); Figure 2i
5b	I	$O_2^{\bullet-}$	Figure 1a (main text); Figure S4
6a	I	$O_2^{\bullet-}$, 1O_2	Figure 4a (main text); Figure 4b,c
6b	I	$O_2^{\bullet-}$, 1O_2	Figure 4a (main text); Figure 4b,c
6c	I	$O_2^{\bullet-}$, 1O_2	Figure 4a (main text); Figure 4b,c


3. AI Image Recognition

Use InDraw's image recognition function to automatically convert structural images into an editable molecular graph.

1) Copy the image to InDraw, then right click and select AI 图像识别 (AI Image Recognition).




2) Click AI 图像识别 (AI Image Recognition) and take a screenshot for recognition.



4. Human-Assisted Inspection

Manually inspect the recognized structure and correct any errors in atoms, bonds, or stereochemistry.



5. Convert Molecules to SMILES

Select the verified molecule, and right click→复制为 (Copy As)→SMILES, or directly use the shortcut keys “Alt+Ctrl+C” to obtain the SMILES string.

